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IMPORTANT : Please read these instructions carefully before use**LIABILITY**

The user of DEEP Company prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any, direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications regularly. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

DESCRIPTION AND INTENDED USE

LIKE prosthetic surfaces are designed to facilitate prosthetic restoration, protect LIKE implants during healing, and ensure the connection between LIKE implants (which replace the root) and the crown. The intended use of each component is given below:

Family product	Intended use
Prosthetic abutment	Connect to the implant in order to support the fixed prosthesis and ensure its junction with the implant
Healing screw	Prepare the tissues during the healing phase
Cover screw	Protect the connectors of the implant during the healing phase

MATERIAL

The prosthetic parts are made of titanium grade 5 (TA6V according to ISO 5832-3):

Element	Mass Fraction %
Aluminum (Al)	5.5 - 6.5
Vanadium (V)	3.5 – 4.5
Iron (Fe)	0.3 max
Oxygen (O)	0.2 max
Nitrogen (N)	0.05 max
Hydrogen (H)	0.015 max
Titanium (Ti)	Balance

Except for the multi-abutment protective cap, which is machined from PEEK:

Element	Mass Fraction %
Polyether ether ketone (PEEK)	100

INDICATIONS

Prosthetic components are used to make prosthetic rehabilitations that range from replacing one single tooth to an entire arch of bridgework, as well as retentive elements for overdenture applications.

INTENDED PATIENT POPULATION

LIKE prosthetic surfaces are intended for individuals missing at least one tooth and who have no contraindications as identified in the following section.

The patient profile is highly varied (in terms of age and gender) because the reasons for implantation can be either traumatic or non-traumatic. Younger patients are more often treated for traumatic reasons, while older patients typically require implants due to diseases affecting the bone and/or periodontal tissues.

CONTRAINDICATIONS

The placement of implants is not recommended in the following cases:

- Allergy or hypersensitivity to the chemical components present in the materials used and mentioned in the "Description and intended use" section ;
- Patient with the following medical risk factors: cardiovascular diseases, diabete mellitus, osteoporosis and other bone diseases requiring bisphosphonates, diseases requiring local or systemic treatments (anticoagulant, chemotherapy, radiotherapy, steroid), severe mental illnesses (e.g., dementia or schizophrenia), psychotic disorders, organ transplant patients, bleeding disorders, an immune deficiency, patients taking some pharmacological agents such as proton pump inhibitors, physical capacity not adapted to implant treatment.
- Patient with the following local contraindications : history of recurrent diseases such as dental and periodontal diseases, unfavorable occlusion for whom it is impossible to place the adequate size or number of implants, bruxism or other parafunctional habits, bad oral condition (presence of cavities or lesions near the tooth loss) and inadequate bone volume, unless an augmentation procedure can be considered
- Patient with the following behaviour : poor oral hygiene, smoking, lack of motivation or cooperation.

In the elderly, implant placement is not contraindicated if they have adequate cognitive and physical abilities. For children, the placement of implants is not recommended until the end of jaw growth has been duly established. Preoperative gum and bone deficiencies can compromise the aesthetic result or result in unfavourable implant positioning. In pregnant women, implant placement is not recommended because it may be associated with taking an x-ray. There is no contraindication for breastfeeding women.

INTENDED USER

LIKE prosthetic surfaces must be handled and/or implanted by dentists, maxillofacial surgeons, or stomatologists trained in specific operating techniques for the placement of dental implants as part of their initial or continuing training, and having read the instructions of use. We strongly recommend that the practitioner always undergo special training before implementing a new treatment method.

The clinician is responsible for any complications that can occurred from a use that do not comply with the indications and contraindications, from a defective surgical technique or a lack of asepsis. These complications shall in no way be attributed to DEEP COMPANY.

EXPECTED CLINICAL BENEFIT

The clinical benefits of the LIKE Implant System product, meaning the implant associated with abutments, are:

- increase of oral health quality of life;
- satisfaction with aesthetic and function.

EXPECTED CLINICAL PERFORMANCE

Survival rate (using Kaplan-Meier survival analysis) is :

- Superior to 98% at 1 year;
- Superior to 95% at 5 years, and;
- Superior to 90% at 10 years;

Expected lifetime of Dental implant system LIKE is 10 years.

COMPLICATIONS AND UNDESIRABLE EFFECTS

- loss of bone,
- loss of implant
- Irritation/inflammation
- infections,
- mechanical failures including fatigue fracture of implants
- modification of the treatment plan
- lengthening of the operating time
- surgical implant explantation
- allergies to different materials
- Healing disorder
- Pain

PATIENT INFORMATION:

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.
- informing the patient that he must report to DEEP COMPANY and the competent health authorities any serious incident occurring with LIKE prosthetic surfaces.

PROCEDURAL PRECAUTIONS

Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw; avoiding excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

INTERACTION WITH MAGNETIC RESONANCE IMAGING AND OTHER ELECTROMAGNETIC SYSTEM

The safety and compatibility of LIKE prosthetic surfaces that remain in the patient's body in an MR environment have not been proven. These devices have not been tested for heat buildup, migration, or artifacts in MR environments.

Performing an MRI scan on a patient with such a device may cause injury, and it is therefore contraindicated.

COMPATIBILITY INFORMATION:

The LIKE A/CC prosthetic surfaces are compatible with LIKE-A and LIKE-CC implant ranges, as well as LIKE-A and LIKE-CC instrument ranges.

The LIKE-BS prosthetic surfaces are compatible with LIKE-BS implants and instruments.

The LIKE-I prosthetic surfaces are compatible with LIKE-I implants and instruments.

PACKAGING, STORAGE AND HANDLING

From delivery until use, store devices in a clean room away from products which could have a corrosive action.

In the event of an apparent packaging defect, the product must be returned to DEEP Company.

PROTOCOL OF CLEANING AND STERILIZATION BEFORE USE

The non-sterile prosthetic surfaces supplied by DEEP COMPANY have been cleaned of manufacturing residues and then washed but do not have a sufficient level of decontamination to allow them to be used directly. Therefore, before first use, all LIKE prosthetic surfaces must be cleaned, disinfected and sterilized according to the protocol mentioned below.

The protocol was drawn up in accordance with the SF2S guide to Good Sterilisation Practice for Reusable Medical Devices, updated in 2021.

Concerned devices	All LIKE prosthetic surfaces are covered by this protocol
Warning	The cleaning, disinfection and sterilisation protocols have been validated by DEEP Company. We therefore recommend that you follow each step scrupulously (parameters, detergent, etc.). For cleaning and disinfection, DEEP Company has used the detergent Neodisher® MediClean Dental (DrWeigert) for the validation of the protocols. However, other mild alkaline protease enzymatic detergent may be used if they are CE marked and aldehyde-free. The use of other detergents is the responsibility of the user and must be validated before use. Protective clothing (mask, gloves and goggles) must be worn when handling cleaning and disinfecting detergent.

Instructions	
1. Automated cleaning & disinfection	Equipment : For this operation, DEEP COMPANY recommends automatic cleaning and disinfection of its dental instruments using a washer-disinfector that complies with standard NF EN ISO 15883-1&-2. The washer/disinfector must be validated beforehand by the user and regularly maintained. Detergent : Neodisher® MediClean Dental detergent (DrWeigert), or mild alkaline protease enzymatic detergent, at a concentration of 0.2%.

DEEP Company thanks you for your support and remains at your disposal for any additional information you may require

DEEP Company
ZA du bois Saint-Pierre
38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03

www.astem.fr

contact@deep-company.com

IMPORTANT : Please read these instructions carefully before use

	Steps : 1. Place devices in a wash basket suitable for the washer-disinfector. Orient in manner that allows water to flow into their cavities. Ensure that all surfaces are exposed to the cleaning solution. ⚠ Do not overload the wash basket. The load is at the user's initiative and must be validated by the user during the initial qualification of the equipment. ⚠ Avoid any contact between the devices, which may cause damage during washing 2. Start the next cycle: <table><tr><th>Steps</th><th>Water*</th><th>Minimal temperature</th><th>Minimal time</th></tr><tr><td>Prewash</td><td>Purified water</td><td>Cold water : < 45 °C</td><td>1 minute</td></tr><tr><td>Washing</td><td rowspan="2">Purified water + detergent</td><td>Heated water : 50-60°C</td><td>10 minutes</td></tr><tr><td>Thermal disinfection</td><td>Heated water : 85°C-90°C</td><td>5 minutes</td></tr><tr><td>Rinsing</td><td>Purified water</td><td>Cold water : <45 °C</td><td>1 minute</td></tr><tr><td>Drying</td><td>N/A</td><td>95°C</td><td>20 minutes</td></tr></table> * Water quality according to European Pharmacopeia in force 3. When unloading, inspect the surfaces and cavities of medical devices to verify complete removal of visible dirt. If necessary, repeat this cleaning cycle. 4. Likewise, check that the prosthetic surfaces do not show any deformation, damage or defects (loss of marking), caused by the cleaning step, and which could affect the safety and performance of the device. If this is the case, please stop the operation and return the defective products to DEEP Company (see "material return" section).	Steps	Water*	Minimal temperature	Minimal time	Prewash	Purified water	Cold water : < 45 °C	1 minute	Washing	Purified water + detergent	Heated water : 50-60°C	10 minutes	Thermal disinfection	Heated water : 85°C-90°C	5 minutes	Rinsing	Purified water	Cold water : <45 °C	1 minute	Drying	N/A	95°C	20 minutes	
Steps	Water*	Minimal temperature	Minimal time																						
Prewash	Purified water	Cold water : < 45 °C	1 minute																						
Washing	Purified water + detergent	Heated water : 50-60°C	10 minutes																						
Thermal disinfection		Heated water : 85°C-90°C	5 minutes																						
Rinsing	Purified water	Cold water : <45 °C	1 minute																						
Drying	N/A	95°C	20 minutes																						
2. Drying	The drying step is integrated into the washing/disinfection step above. Visually check that the devices are dry when unloading the devices. If this is not the case, carefully wipe all device surfaces with a clean, lint-free cloth (single use) until the prosthetic surface are dry.																								
3. Packaging	Type of packaging to be used : DEEP COMPANY recommends the use of a pre-formed sterile barrier system, composed of a sealed external and internal sterilisation sheaths, that complies with standard NF EN ISO 11 607-1 and is suitable for autoclave sterilisation. In addition, these sheaths must be adapted to the dimensions of the surgical kit. ⚠ Please check the condition of the sterilisation sheaths before packaging. Steps : 1. Place the prosthetic surfaces in the appropriate sterilization pouch. 2. Check that the packaging is properly closed, then affix a traceability label containing the following information: the product reference and IUD, the type of sterilisation to be used, the use-by date (depending on the sachet used) and the applicable storage conditions if applicable.																								
4. Moist Heat Sterilization	Equipment: DEEP COMPANY recommends performing the sterilisation of its dental instruments using a CE-marked, validated autoclave that complies with standard NF EN ISO 17665-1. The autoclave must be regularly maintained in order to maintain sterilization efficiency. Steps : 1. Place the surgical kit, previously packaged, in the autoclave basket to ensure even steam penetration through the packaging. Make sure there is no contact between the packaging and the wall of the steriliser; ⚠ Do not overload the autoclave basket. The load is the user's responsibility and must be validated by the user during the initial qualification of the equipment. 2. To check the effectiveness of the sterilisation cycle, place a physical-chemical indicator on the load to be sterilised. 3. Start a sterilisation cycle using the parameters below *: <table><tr><th>Parameters</th><th colspan="2">Values</th></tr><tr><td>Cycle type</td><td colspan="2">Prevacuum</td></tr><tr><td rowspan="4">Pre-treatment</td><td colspan="2">3 negative pulses (< 100 mbar)</td></tr><tr><td colspan="2">2 positive pulses (1500 ± 100 mbar)</td></tr><tr><td colspan="2">2 positive pulses (1900 ± 100 mbar)</td></tr><tr><td colspan="2">2 negative pulses (1000 ± 100 mbar)</td></tr><tr><td>Set point temperature</td><td>France 134°C</td><td>Europe 134°C</td></tr><tr><td>Exposure time</td><td>18 min</td><td>3 min and 8 sec</td></tr><tr><td>Dry time under vacuum</td><td colspan="2">30 min</td></tr></table> * The choice of cycle must comply with the standards and regulations in force in the country of sale 4. After unloading the load, check the physical-chemical indicator (colour change of sterilisation indicators) and the sterilisation efficiency (integrity of the packaging, no humidity, integrity of the traceability label, digital recordings of the	Parameters	Values		Cycle type	Prevacuum		Pre-treatment	3 negative pulses (< 100 mbar)		2 positive pulses (1500 ± 100 mbar)		2 positive pulses (1900 ± 100 mbar)		2 negative pulses (1000 ± 100 mbar)		Set point temperature	France 134°C	Europe 134°C	Exposure time	18 min	3 min and 8 sec	Dry time under vacuum	30 min	
Parameters	Values																								
Cycle type	Prevacuum																								
Pre-treatment	3 negative pulses (< 100 mbar)																								
	2 positive pulses (1500 ± 100 mbar)																								
	2 positive pulses (1900 ± 100 mbar)																								
	2 negative pulses (1000 ± 100 mbar)																								
Set point temperature	France 134°C	Europe 134°C																							
Exposure time	18 min	3 min and 8 sec																							
Dry time under vacuum	30 min																								

	cycle parameters). Keep the indicator and the checks in the patient file to ensure traceability.
5. Storage	Please keep the prosthetic surfaces in their sterilization packaging in a dry and clean environment until the expiration date (inscribed on the traceability label). Please consult the notice of packaging used for further information on storage conditions, if applicable. Before use, please check the packaging and instruments (integrity of the packaging, sign of humidity, etc.). In the event of damage to the packaging, sterility cannot be guaranteed.

DEVICE REPROCESSING

LIKE prosthetic surfaces are single use devices. Reuse could lead to crossed contamination.

USE:

The application and removal of DEEP Company products must be carried out using appropriate and compatible instruments (see chapter COMPATIBILITY INFORMATION).

Guide to tightening torques:

Surface prosthetic	Tightening torques (N.cm)	
	LIKE-I	LIKE-BS/-CC / -A
Direct implant monobloc abutments (without removable screw)	20	35
Direct implant screw-retained abutments	15	35
Cover and healing screws	Manual	
Indirect screw-retained abutments (prosthetic surfaces going on the multi abutments)	15	

PRODUCT DISPOSAL

Sorting between infectious waste and waste comparable to household waste lies with the professional who produced the waste.

- Typically common waste products such as packaging are similar to household waste and follow the same waste management system.
- Items and materials soiled by blood and / or body fluids have to be treated as waste from risky treatment activities.

PRODUCT RETURN

In the event of device return following an order error, it must be returned unused in its original unopened packaging.

All incidents related to the use of LIKE prosthetic components must imperatively be reported to DEEP Company.

Moreo

Additionally, any serious incident related to the use of LIKE instruments must be reported to the relevant health authorities.

Please clean and sterilize the device(s) according to the previous described procedures and send them in double sterile pouch. The return will be accompanied by the " Implant removal management" form available in the "documentary" section of our website

ADDITIONNAL INFORMATION

For more information on the indication and use of different systems, please contact DEEP Company by the following email : contact@deep-company.com or on our website : www.astem.fr

Summary of the safety and performance characteristics of LIKE Prosthetic Components are available in the EUDAMED database, accessible at the following address: <https://ec.europa.eu/tools/eudamed>.

While waiting for EUDAMED to become operational, summary of the safety and performance characteristics for LIKE prosthetic surfaces can be requested via email at contact@deep-company.com.

MEANING OF THE SYMBOLS

The meaning of the symbols available on the external label is available in the following table :

	Legal manufacturer		Quantity
	Material		Reference
	Medical device		Do not reuse
	Manufacturing date		Unique Device Identifier
	Batch number		Do not use if packaging is damaged
	Consult electronic instructions for use		CE marked

DEEP Company thanks you for your support and remains at your disposal for any additional information you may require

DEEP Company
ZA du bois Saint-Pierre
38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03

www.astem.fr

contact@deep-company.com