### RESPONSABILITES :

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### LIABILITY

The user of LIKE implants has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP COMPANY disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP COMPANY products. The clinician is also obliged to study the latest developments in regard to the DEEP COMPANY products and their applications. In cases of doubt, the clinician has to contact DEEP COMPANY. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP COMPANY assumes no liability whatsoever for damages arising thereof.

### DESCRIPTION

Endosseous threaded sandblasted dental implant.

### INTENDED PURPOSE

LIKE implants are intended to replace the root of an absent tooth, in the mandible as in the maxilla, caused by accidental or pathological loss, extraction or agenesis, when the other medical and surgical possibilities have been considered less appropriate.

### MATERIAL

LIKE implants are made from pure titanium grade 4B according to ISO 5832-2 standard :

Element	Mass fraction %
Nitrogen (N)	0.05 max
Carbon (C)	0.08 max
Hydrogen (H)	0.012 max
Iron (Fe)	0.5
Oxygen (O)	0.4
Titanium (Ti)	The rest

### INDICATIONS

Dental implants are used for replacing teeth both in mandible and maxilla. They must be placed in the jaw bone to be integrated (osteointegration). Restorations range from a single replacement to the use of elements for overdenture applications as well as the full arch of bridgework.

### INTENDED PATIENT POPULATION

LIKE implants are intended for people with the absence of at least one tooth who does not have any contraindications identified in the instructions for use.

The patient population targeted by these devices is very varied (age, sex) because the causes leading to implantation can be traumatic or atraumatic. Young patients will be treated more for traumatic causes and elderly patients more for pathologies affecting bone and / or periodontal tissues.

### CONTRAINDICATIONS

The placement of implants is not recommended in the following cases:

- Allergy or hypersensitivity to the chemical components present in the materials used and mentioned in the "Description and intended use" section ;

- Patient with the following medical risk factors: cardiovascular diseases, diabete mellitus, osteoporosis and other bone diseases requiring bisphosphonates, diseases requiring local or systemic treatments (anticoagulant, chemotherapy, radiotherapy, steroid), severe mental illnesses (e.g., dementia or schizophrenia), psychotic disorders, organ transplant patients, bleeding disorders, an immune deficiency, patients taking some pharmacological agents such as proton pump inhibitors, physical capacity not adapted to implant treatment.

- Patient with the following local contraindications : history of recurrent diseases such as dental and periodontal diseases, unfavorable occlusion for whom it is impossible to place the adequate size or number of implants, bruxism or other parafunctional habits, bad oral condition (presence of cavities or lesions near the tooth loss) and inadequate bone volume, unless an augmentation procedure can be considered

- Patient with the following behaviour : poor oral hygiene, smoking, drug/alcohol abuse, lack of motivation or cooperation.

In the elderly, implant placement is not contraindicated if they have adequate cognitive and physical abilities. For children, the placement of implants is not recommended until the end of jaw growth has been duly established. Preoperative gum and bone deficiencies can compromise the aesthetic result or result in unfavourable implant positioning. In pregnant women, implant placement is not recommended because it may be associated with taking an x-ray. There is no contraindication for breastfeeding women.

### INTENDED USERS

LIKE implants must be handled and/or implanted by dentists, maxillofacial surgeons, or stomatologists trained in specific operating techniques for the placement of dental implants as part of their initial or continuing training, and having read the instructions of use. We strongly recommend that the practitioner always undergo special training before implementing a new treatment method.

The clinician is responsible for any complications that can occur from a use that do not comply with the indications and contraindications, from a defective surgical technique or a lack of asepsis. These complications shall in no way be attributed to DEEP COMPANY.

### EXPECTED CLINICAL BENEFIT

The clinical benefits of the LIKE Implant System product, meaning the implant associated with abutments, are:

- increase of oral health quality of life;
- satisfaction with aesthetic and function.

### EXPECTED CLINICAL PERFORMANCE

Survival rate (using Kaplan-Meier survival analysis) is :

- Superior to 98% at 1 year;
- Superior to 95% at 5 years; and;
- Superior to 90% at 10 years;

Expected lifetime of Dental implant system LIKE is 10 years.

### COMPLICATIONS AND UNDESIRABLE EFFECTS

- loss of bone,
- loss of implant
- Irritation/inflammation
- infections,
- mechanical failures including fatigue fracture of implants
- modification of the treatment plan
- lengthening of the operating time
- surgical implant explantation
- allergies to different materials
- Healing disorder
- Pain

### PRE-OPERATIVE PRECAUTIONS

A careful clinical and radiological examination of the patient has to be performed prior to surgery to determine the physical and psychological status of the patient. Special attention has to be given to patients who have local or systemic factors that could interfere with the healing process of either bone or soft tissue or the osteointegration process (cigarette smoking, uncontrolled diabetes, oro-facial radiotherapy, steroid therapy, infections in the neighboring bone). In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism or unfavorable jaw relationships, reappraisal of the treatment option may be considered.

### PER-OPERATIVE PRECAUTIONS

Particular caution should be used when placing small diameter implant (LIKE-A, LIKE-I, LIKE-CC : Ø3.5, LIKE-BS : Ø3.3 NP) in the posterior region due to the risk of prosthetic overload. Lack of adequate quantity and/or quality of remaining bone, infection and generalized diseases may be potential causes for failure of osteointegration both immediately or after the surgery. Besides the mandatory precautions for any surgery (such as asepsis), during drilling in the jaw bone, the clinician must avoid damage the nerves and vessels by referring to anatomical knowledge and preoperative radiographs. Failure to recognize actual lengths of drills relative to radiographic measurement or drilling beyond the depth intended can result in permanent injury to nerves or other vital structures. It can result in a permanent numbness to the lower lip and chin or lead to hemorrhage such as in the floor of the mouth.

Implants may be tilted up to 45° relative to the occlusal plane. When used with angulations between 30° and 45°, the following applies: The tilted implants must be splinted; a minimum of 4 implants must be used when supporting a fixed prosthesis in a fully edentulous arch.

The practitioner is advised to provide regular patient follow up and to inform about good oral hygiene.

### PATIENT INFORMATION

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the safety and effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.
- informing the patient that he must report to DEEP COMPANY and the competent health authorities any serious incident occurring with LIKE implants.

### INTERACTION WITH MAGNETIC RESONANCE IMAGING AND OTHER ELECTROMAGNETIC SYSTEM

The raised temperature after interaction between MRI and implants does not pose any known hazard to the patient.

But the quality of the MR image may be compromised, if the area of interest is in the same area of the device: an artefact may appear on the picture.

### STERILITY

All implants are delivered sterile (by gamma ray) and must be used prior to the labelled expiration date.

Warning : Do not use implants if the packaging has been damaged or previously opened.

They are for a single use only. Never reprocess an implant. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination. In case of re-sterilization of the product, DEEP COMPANY declines all responsibility regardless of the person who carried out the re-sterilization or the method used.

### SURGICAL PROCEDURES

- During drilling procedures, bone quality should be considered (see tables below). Drilling must proceed at high speed (maximum 1200 rpm for trocar drill and 800 rpm for step/twist drills) under constant and profuse irrigation by sterile saline at room temperature.
- Prepare implant site. When using a flapless surgery, add-on soft tissue height to drill depth.
- Measure the final depth of implant site for applicable implant length using a depth probe. Use the drills graduations to perform a drilling with the proper depth.  
**Caution:** Take into account the additional drill length when drilling near vital structures.
- Open the implant package and pick up the implant from inner casing with the implant driver. The implants are ideally installed with low speed, max 25 rpm, using a surgical machine or a manual torque wrench.
- Place and tighten the implant without exceeding an insertion torque of 45 N.cm (LIKE-BS, LIKE-CC, LIKE-I) or 70 N.cm (LIKE-A) maximum.  
**Caution:** Never exceed the insertion torque of 45 N.cm (LIKE-BS, LIKE-CC, LIKE-I) or 70 N.cm (LIKE-A). Overtightening an implant may lead to damage of the implant, fracture or necrosis of the bone site.
- If the implant gets stuck during implant installation or 45 N.cm (LIKE-BS, LIKE-CC, LIKE-I) or 70 N.cm (LIKE-A) is achieved before fully seated, rotate the implant counter clockwise using drilling machine or manual torque wrench, in reverse mode, and remove the implant from the site. Replace the implant back into titanium casing. Use a wider drill, countersink or screwtap to enlarge the site. If screwtap is used, place it into the implantation site, select the low speed mode (25 rpm) and drill to the appropriate length. Switch to reverse mode and pull the screwtap out. Continue the implant installation until desired position is achieved.
- For the immediate function, the implant should be able to withstand a minimal final torque of 35 Ncm.
- Depending on surgical protocol, place a cover screw or an abutment and suture.

For any additional information on surgical procedures please consult the "LIKE product documentations" treatment guidelines available at <https://astem.fr/documentation/>.

### Recommended drill sequences based on bone quality.

Note: all data is stated in mm.

Drills within brackets (--) denote widening of the cortical bone only, not drilling the full drilling depth.

Table A : LIKE-A drill sequence

LIKE-A Implant diameter	soft bone Type IV	medium bone Type II-III	dense bone Type I
ø 3.5	2.0 (2.4/2.8)	2.0 2.4/2.8 (2.8/3.2)	2.0 2.4/2.8 (2.8/3.2)
ø 4.3	2.0 2.4/2.8 (2.8/3.2)	2.0 2.4/2.8 (3.2/3.6)	2.0 2.4/2.8 (3.2/3.6) (3.8/4.2)
ø 5.0	2.0 2.4/2.8 3.2/3.6	2.0 2.4/2.8 3.2/3.6 3.8/4.2	2.0 2.4/2.8 3.2/3.6 3.8/4.2 (4.2/4.6)

Table B : LIKE-BS drill sequence

LIKE-BS implant diameter	Soft bone Type IV	medium bone Type II-III	dense bone Type I
ø 3.3	2.0	2.0	2.0 2.4/2.8
ø 3.75	2.0 (2.4/2.8)	2.0 2.4/2.8 3.0	2.0 2.4/2.8 3.2
ø 4.0	2.0 (2.4/2.8)	2.0 2.4/2.8 3.2	2.0 2.4/2.8 3.4
ø 5.0	2.0 2.4/2.8 3.0	2.0 2.4/2.8 3.2/3.6	2.0 2.4/2.8 3.2/3.6 3.8/4.2

Table C : LIKE-CC drill sequence

LIKE-CC implant diameter	Soft bone Type IV	medium bone Type II-III	dense bone Type I
ø 3.5 NP	Trocar 2.0 3.5	Trocar 2.0 3.5	Trocar 2.0 3.5 (4.3)
ø 4.3 RP	Trocar 2.0 3.5 4.3	Trocar 2.0 3.5 4.3	Trocar 2.0 3.5 4.3 (5.0)
ø 5.0 WP	Trocar 2.0 3.5 4.3 5.0	Trocar 2.0 3.5 4.3 5.0	Trocar 2.0 3.5 4.3 5.0

Table D : LIKE-I drill sequence

LIKE-I implant diameter	Soft bone Type IV	medium bone Type II-III	dense bone Type I
ø 3.5	Trocar 2 (2/3)	Trocar 2 2/3	Trocar 2 2/3 Terminal drill 3.53.5
ø 4	Trocar 2 2/3	Trocar 2 2/3 Terminal drill 3.5	Trocar 2 2/3 Terminal drill 3.54
ø 4.5	Trocar 2 2/3 Terminal drill 3.5	Trocar 2 2/3 3/4	Trocar 2 2/3 3/4 Terminal drill 3.54.5

### INFORMATION ON COMPATIBILITY

LIKE A-implants are compatible with LIKE-A prosthetic parts and instruments

LIKE-BS implants are compatible with LIKE-BS prosthetic parts and instruments

LIKE-CC implants are compatible with LIKE-CC prosthetic parts and instruments

LIKE-I implants are compatible with LIKE-I prosthetic parts and instruments

### ADVERSE EVENT AND RETURN

All incidents related to the use of LIKE implants must imperatively be reported to DEEP COMPANY. Moreover, any serious incident linked to the use of DESIGN implants must be reported to the competent health authorities.

Please clean and sterilize the device(s) according to the procedure described in the instruction for use of LIKE instruments or prosthetics parts and return the devices under double pouch. The return will be accompanied by the "information collection" form available in the "documentation" section of our website : <https://astem.fr/documentation/>

### ADDITIONNAL INFORMATION

For more information on the indication and use of different systems, please contact DEEP COMPANY by the following email : [contact@deep-company.com](mailto:contact@deep-company.com) or on our website : [www.astemdigital.fr](http://www.astemdigital.fr)

Summaries of the safety and performance characteristics of LIKE implants are available in the EUDAMED database, accessible at the following address: <https://ec.europa.eu/tools/eudamed>. Until EUDAMED is operational, the summary of performance and safety of LIKE implants are available on our website, in the "documentation" section of our website : <https://astem.fr/documentation/>.

### MEANING OF THE SYMBOLS :

The meaning of the symbols available on blister and external labels are reported in the following table :

	Manufacturer		Unique sterile barrier system with inner protective packaging
	Expiry date		Reference
	Sterilized by irradiation		Do not reuse
	Manufacturing date		Unique device identifier
	Lot		Do not use if packaging is damaged
	Do not resterilize		CE marked by TÜV SÜD
	Medical device		Consult electronic instructions for use